

Headspace Gas Chromatography (HS-GC)

How Safe is Your Mask?




Ethylene Oxide

Surgical masks are typically sterilized with ethylene oxide at room temperature. However, ethylene oxide is a carcinogen. There are strict requirements on the ethylene oxide residue in single-use sterile medical device products. Both ISO 10993.7-2008 (Biological evaluation of medical devices) and China's National standard "Biological Evaluation of Medical Devices Part VII: Ethylene Oxide Sterilization Residues" (GB / T16886.7-2015) specify the maximum allowable residual amount of ethylene oxide. Shimadzu HS-GC is most suited for measurement of ethylene oxide in safety masks

Easy to use | Minimal sample preparation | High performance | Rugged

Headspace Gas-Chromatography (HS-GC)

Methodology involves simple sample preparation of transferring 1 g of sample into headspace bottle followed by sealing for testing. Excellent linearity was achieved for calibration curve ranging from 0.2-10 ug standard, along with reproducibility of 2.64% for 6 injections of 0.4 ug ethylene oxide.



Chromatogram showing ethylene oxide measurement in standard sample, with calibration curve of 0.2 – 10 ug in the insert.

GC Condition		HS Condition	
Column	SH-Rxi™-624Sil MS 60 m, 0.32 mm ID, df=1.8 µm (P/N: 221-75963-60)	Equilb. temp.	70 °C
Oven temp. program	40°C (3 min) @ 5 °C/min to 85°C @ 30 °C/min to 200°C (2 min)	Sample line temp.	100 °C
Injector	Split @ 220 °C ; Ratio 10:1	Transfer line temp.	110 °C
Detector	TCD @ 250 °C	Equilb. time	30 min
		Load time	0.5 min

HS-10 with GC-2014 for versatile configuration



※For liquid injection, please remove HS-10 transfer line and set the AOC-20i

Switch between headspace and liquid injections with HS-Switcher for high accuracy and precision in various applications.

HS-20 with Nexis GC-2030 for highest throughput



With the shortest inert transfer line in the industry (30 cm), we ensure minimal carryover between different analyses.

The headspace autosampler also recognizes 10 mL and 20 mL automatically. It hence allows you to prepare various samples without the need for adapters.

Versatility with varied samples	Great for	High-throughput routine labs
20 vials	Headspace Vials	90 vials
6 vials	Simultaneous Heating	12 vials
225 °C	Oven Temperature	300 °C
225 °C	Loop Temperature	300 °C
225 °C	Transfer Line Temperature	350 °C
Yes	Vial Shaker	Yes
No	Leak Check	Yes
No	Barcode Reader	Yes

Shimadzu GC column (SH-Rxi™-624Sil MS)

ID	df	Temp. Range	20m	30m	60m	75m	105m
0.18 mm	1.00 µm	-20 to 300/320 °C	227-36075-01	–	–	–	–
0.25 mm	1.40 µm	-20 to 300/320 °C	–	221-75962-30	227-36076-01	–	–
0.32 mm	1.80 µm	-20 to 300/320 °C	–	227-36077-01	221-75963-60	–	–
0.53 mm	3.00 µm	-20 to 280/300 °C	–	227-36078-01	227-36078-02	227-36078-03	227-36078-04



Shimadzu Europa
www.shimadzu.eu

For Research Use Only. Not for use in diagnostic procedures.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

Company names, product/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation or its affiliates, whether or not they are used with trademark symbol "TM" or "®". Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services. Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The contents of this publication are provided to you "as is" without warranty of any kind, and are subject to change without notice. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication.