

Technical Report

RF-20A Fluorescence Detector Basics and Applications

Abstract:

A fluorescence detector is one of the typical detectors in an HPLC analysis. Its operational principle and characteristics enable a highly sensitive and selective analysis compared to a UV absorbance detector which is the first choice of detector in HPLC. Although the number of compounds that can be fluorescently detected is limited, a derivatization makes the technique applicable to a wide range of analyses. This leads to a highly sensitive quantitative analysis in many fields such as pharmaceutical, food, and environmental markets. In this report, the principle of the fluorescence detector and its advantages are described. Additionally, example analyses of Shimadzu fluorescence detector RF-20Axs will be introduced to showcase its unique temperature control function (e.g., cooling).

1. Basics of Fluorescence Detection Technique

When light of a specific wavelength (excitation wavelength; Ex.) is irradiated onto a substance, that substance emits light that usually has a longer wavelength (emission wavelength; Em.), and is referred to as fluorescence. Unlike an ultraviolet (UV) detector, which measures the amount of light absorbed by a substance at a specific wavelength, the light that is emitted from the sample due to fluorescence is split into a spectrum by a monochromator, and the intensity of light at a specific emission wavelength is measured. Fluorescence detection has the following features:

- Selectivity is high because measurement is conducted using specific excitation and emission wavelengths specific to target substance.
- Fluorescent substances can be detected with high sensitivity.
- 3) Sensitivity and selectivity can be improved by using the derivatization techniques.

Fluorescence detection is generally used for analysis when sensitivity and selectivity are required, especially when the analyte has little or no UV absorbance and can be derivatized to produce fluorescence. Although LC/MS is increasingly being used for high-sensitivity analysis, fluorescence detection by HPLC is often the officially regulated method. In addition, compared to LC/MS, fluorescence detection is more economical and maintenance is easier. Here we describe the fundamental principles of the fluorescence detector, and present actual examples of analysis using the Shimadzu fluorescence detector RF-20Axs to illustrate the innovations that provide stable and high sensitivity analysis.

1-1. Principle of Fluorescence

Fluorescent substances conform to the following principle of fluorescence. The substance enters a state of excitation when it is irradiated with excitation light. Due to the instability of this excited state, the substance quickly returns to its original state (ground state). The energy emitted at this time is referred to as fluorescence, and a fluorescence detector measures the intensity of the emitted fluorescence. Since excitation is conducted at a specific wavelength, and measurement of the fluorescence is conducted at a specific wavelength, selectivity is higher than with absorbance and other types of detectors (Fig. 1).





>>>> Learn More

Analysis Basics

• Basics of LC

• Fluorescence Detection

Fig. 1 Mechanism of Fluorescence from a Substance

1-2. Characteristics of Fluorescence Emission

Fluorescence generates a spectrum like that shown in Fig. 2. It is important to note that upon irradiation with the excitation wavelength, excitation scattered light and Raman scattered light will always occur in the vicinity of the excitation wavelength. Therefore, if the excitation wavelength and fluorescence (emission) wavelength are too close to one another, it may be difficult to obtain a valid chromatogram. In addition, background noise will increase if fluorescence is generated by the solvent. For this reason, it is important to use HPLC-grade water and solvent for the mobile phase.



Fig. 2 Characteristics of Fluorescence Spectrum

Detection Method	Typical Sensitivity
Absorbance	ng
Fluorescence	pg
Differential refractive index	μg
Electrochemical	pg
Electroconductivity	ng

Table 1 Typical Levels of Sensitivity of Various Detection Methods



Shimadzu Fluorescence Detector: RF-20Axs

1-3. Sensitivity and Selectivity of Fluorescence Detectors

As explained in the principle of fluorescence, two settings are made in the fluorescence detector : the excitation wavelength and the emission wavelength. This is what allows the fluorescence detector to provide greater selectivity than other detectors. The data below are chromatograms obtained from analysis of vitamin B₂ in a food product. In the chromatogram of Fig. 3 obtained using absorbance detection, the contaminant peaks overlap the analyte peak due to low selectivity. Furthermore, the weak UV absorbance intensity of vitamin B₂ makes quantitation difficult. On the other hand, as shown in Fig. 4, a fluorescence detector provides highly selective detection in vitamin B₂ analysis, because contaminants in the sample are not detected at the excitation and fluorescence wavelengths used for vitamin B_2 detection. In addition, since the background can be suppressed to a low level by using a non-fluorescent mobile phase, detection with high selectivity as well as high sensitivity is possible. The typical levels of sensitivity of various detection methods are shown in Table 1.



Fig. 3 Absorbance Detection Technique



Fig. 4 Fluorescence Detection Technique

2. Factors that Influence Fluorescence Detection

2-1. Effect of Dissolved Oxygen in Mobile Phase

When the mobile phase contains dissolved oxygen, fluorescent intensity decreases. This is referred to as "quenching." In particular, the intensity may change due to unstable degassing of the mobile phase or prolonged analysis using only offline degassing, either of which could adversely affect peak area repeatability (Fig. 5). Stable degassing can be maintained using the DGU-20A3/As on-line degassing unit in addition to offline ultrasonic degassing to allow analyses that are highly repeatable with excellent sensitivity.



Fig. 5 Change in Peak Intensity due to Dissolved Oxygen in Mobile Phase

2-2. Examples of Temperature Influence in Amino Acid Analysis

We conducted actual amino acid analysis using detector cell temperatures of 25°C, 30°C, 35°C, and 40°C. The chromatogram below (Fig. 6) shows the overlaid results at each temperature. Fig. 7 plots the relative area values (%) for each temperature assuming area values of 100% at 25°C. The peak area values clearly decrease as the detector temperature increases. By increasing cell temperature from 25°C to 40°C, the fluorescent intensities of some of the amino acids decrease by 25% or more. Even when comparing the results at 25°C and 30°C, which is closer to an actual analysis environment, there are substances like Ala which show a decrease in area value of 10% or more, emphasizing the importance of maintaining a constant detector cell temperature from the standpoints of both sensitivity and repeatability.



Fig. 6 Cell Temperature Influence in Amino Acid Analysis



Fig. 7 Relative Area Values (%) at Each Temperature Assuming 100% at 25°C



Nexera[™] Amino Acid Analysis System

3. Features of RF-20Axs Fluorescence Detector

Shimadzu fluorescence detector: RF-20Axs is equipped with a temperature-controlled cell to ensure stable analysis even if the ambient temperature fluctuates. Further, temperature control (with a cooling function) is provided not only for the detector cell, but for the photomultiplier tube (PMT) as well.

Thus, even if the ambient temperature fluctuates widely, the detector is maintained at a constant, near-ambient temperature which provides excellent repeatability without loss of sensitivity. In addition to the improved sensitivity and stability provided with this temperature control in the vicinity of ambient temperature, even further improved sensitivity is achieved with the newly designed optical system and detector cell in the RF-20Axs. As a result, the quantity of light through the flow cell is about 100 times that compared with the previous RF-10AxL model, offering a water Raman *S/N* ratio of at least 2000 for the RF-20Axs. These features are what make this detector a powerful tool, especially when trace-level detection is required.

3-1. Ultra-High Sensitivity Analysis of Anthracene

Fig. 8 shows the results of trace-level analysis of anthracene. An *S/N* ratio of 21.5 was achieved with an injection of 10.48 fg anthracene (RF-20Axs). This is equivalent to a limit of detection of about 1.5 fg (*S/N* ratio = 3), demonstrating excellent sensitivity.

Table 2 Analytical Conditions for Anthracene Analysis

Mobile Phase	Water / acetonitrile = 30 / 70 (v/v)
Flow Rate	0.6 mL/min
Column	Shim-pack™ XR-ODS (50 mm × 2.0 mm I.D., 2.2 µm)
Temperature	40°C
Detection	Ex.: 250 nm, Em.: 400 nm



Fig. 8 Ultra-High Sensitivity Analysis of Anthracene Using RF-20Axs

3-2. Temperature Control (with Cooling) for Highly Stable Analysis

Fig. 9 shows chromatograms that illustrate the effectiveness of cell temperature control when the ambient temperature is fluctuating. By controlling the temperature of the cell, peak intensity barely changes even when the ambient temperature fluctuates. Table 4 shows just how much change occurs in the peak area value at 30°C after the temperature is increased from 25°C relative to a standard area value at 25°C. The peak area repeatability values when the temperature is increased from 25°C to 30°C are also shown. As indicated in the results, excellent repeatability is achieved with the RF-20Axs and its cell temperature.

Table 3 Analytical Conditions for Acridine Analysis

Mobile Phase	Water / acetonitrile = $30 / 70 (v/v)$
Flow Rate	0.6 mL/min
Column	Shim-pack XR-ODS (50 mm × 2.0 mm I.D., 2.2 µm)
Temperature	40°C
Detection	Ex.: 250 nm, Em.: 400 nm
Sample	Acridine



Fig. 9 Comparison of Chromatograms With and Without Cell Temperature Control

 Table 4
 Effectiveness of Cell Temperature Control

 During Ambient Temperature Fluctuation

	Peak Area Change (%)	%RSD During Temperature Fluctuation	
With Cell Temperature Control	0.64	0.29	
Without Cell Temperature Control	-17.45	6.30	

4. Support for Ultra Fast Analysis

Fast response is required to keep up with the sharp peaks obtained in ultra fast analysis. The 10 ms (100 Hz) response of the RF-20A/20Axs permits ultra fast analysis without any loss of resolution. Also, to prevent diminished peak separation due to peak diffusion outside the column during ultra fast analysis, a detector cell having a small capacity is used. The user must be aware of the fact that in fluorescence detection, sensitivity depends on the capacity of the cell. The ultra-high sensitivity capability built into the RF-20Axs permits high sensitivity analysis even during ultra fast analysis with a semi-micro cell. With the RF-20Axs, one unit is all that is required to obtain high performance in all types of analysis, from conventional to ultra fast analysis from conventional to ultra fast analysis.





Table 5 Analytical Conditions for Conventional LC

Mobile Phase	Hexane / 2-propanol = 100 / 0.5 (v/v)
Flow Rate	1.0 mL/min
Column	Shim-pack CLC-SIL (M) (150 mm × 4.6 mm l.D., 5 μm)
Temperature	30°C
Detection	Ex.: 298 nm, Em.: 325 nm

Table 6 Analytical Conditions for Ultra Fast LC

Mobile Phase	Hexane / 2-propanol = $100 / 0.5 (v/v)$
Flow Rate	1.2 mL/min
Column	Shim-pack XR-SIL (75 mmL. × 3 mm I.D., 2.2 µm)
Temperature	30°C
Detection	Ex.: 298 nm, Em.: 325 nm

5. Application Examples

5-1. High-Sensitivity Analysis of Enrofloxacin

Enrofloxacin is one of the new quinolone synthetic anti-bacterial agents (Fig. 11). As a veterinary pharmaceutical product, it is used to treat bovine and swine pneumonia, and as a preventative and remedy for *E. coli* diarrheal syndrome. In Japan, it is not permitted for use with farm-raised fish, and all lots of specific imported products are designated for inspection. Here we conducted high-sensitivity analysis of enrofloxacin using the officially specified analytical method, and we also measured the limit of detection (LOD). The results of analysis using the RF-20Axs indicate an LOD of 0.05 µg/L, as shown in Fig. 12, demonstrating that high-sensitivity detection is possible. In addition, Fig. 13 shows the calibration curve generated using concentrations from 0.05 µg/L to 1000 µg/L. This clearly shows the wide dynamic range, achieving high repeatability and linearity from ultra-trace levels to high concentrations.



Fig. 11 Structural Formula of Enrofloxacin



Fig. 12 Chromatograms of Enrofloxacin at Different Concentration

Table / Analytical Conditions for Enrotioxacin Analy	able /	e / Analvtical	Conditions	tor	Enrofloxacin	i Analvs
--	--------	----------------	------------	-----	--------------	----------

Mobile Phase	Mcilvaine buffer solution (pH 3.0) / acetonitrile = $85 / 15 (v/v)$
Flow Rate	0.5 mL/min
Column	Shim-pack XR-ODS (75 mm × 2.0 mm I.D., 2.2 µm)
Temperature	40°C
Detection	Ex.: 285 nm, Em.: 460 nm
Injection Vol.	2 μL



(0.05 to 1000 mg/L)



5-2. Conventional and Ultra Fast Analysis of Anionic Surfactant (LAS)

We added 0.04 mg/L of LAS to tap water, and analyzed the sample as is without concentrating the sample. The Japanese Ministerial Ordinance for Water Quality Standard specifies a criterion value for a total of 5 LAS compounds (C_{10} to C_{14} in Fig. 14, each 0.04 mg/L, for a total of 0.2 mg/L). The test method specifies HPLC for analysis of the water sample, using a 250× concentration of the solution. However, in this analysis, it is clear that a direct injection can be used without performing any kind of concentration operation. In addition, the time savings achieved using ultra fast analysis allowed the analysis time to be shortened to about 1/5.



Fig. 14 Comparison of LAS Analysis by Conventional Analysis and Ultra Fast Analysis

Table 8 Analytical Conditions for Conventional LC analysis

Mobile Phase	Water / acetonitrile = 35 / 65 (v/v) (including 0.1 mol/L sodium perchlorate)
Flow Rate	1.0 mL/min
Column	Shim-pack VP-ODS (250 mm × 4.6 mm I.D., 4.6 µm)
Temperature	40°C
Detection	Ex.: 221 nm, Em.: 284 nm
Sample	Anionic surfactant standard solution, 20 µL injected

Table 9 Analytical Conditions for Ultra Fast LC Analysis

Mobile Phase	Water / acetonitrile = 35 / 65 (v/v) (including 0.1 mol/L sodium perchlorate)
Flow Rate	1.0 mL/min
Column	Shim-pack VP-ODS (100 mm × 3.0 mm I.D., 4.6 µm)
Temperature	40°C
Detection	Ex.: 221 nm, Em.: 284 nm
Sample	Anionic surfactant standard solution, 8 µL injected

>>>> Learn More

- Analysis of Anionic Surfactants by Prominence[™]-i and RF-20Axs Fluorescence Detector
- High Speed Analysis of Anionic Surfactants

5-3. Simultaneous Analysis of Polycyclic Aromatic Hydrocarbons (PAHs)

The United States Environmental Protection Agency (EPA) specifies 16 polycyclic aromatic hydrocarbons (PAHs) as "priority pollutants," and designates them as restricted substances. Here we conducted multi-wavelength, simultaneous analysis of 15 of these 16 substances. The optimum excitation and emission wavelengths are different for each of the PAH substances, making it necessary to conduct wavelength switching during analysis. However, in high-speed analysis using multiple wavelengths, on-time wavelength switching becomes extremely critical due to the narrow intervals between peaks. The 100 Hz high-speed sampling provided with the RF-20A/20Axs detectors reliably ensures on-time switching, demonstrating their powerful effectiveness in multi-wavelength, high-speed analysis of multiple compounds such as these.



Fig. 15 Simultaneous High-Speed Analysis of 15 Polycyclic Aromatic Hydrocarbons

Table 10 Analytical Conditions for Polycyclic Aromatic Hydrocarbons

Mobile Phase	A: Water, B: acetonitrile
Gradient	B Conc. = 50% (0–0.5 min), 50–88% (0.5–3.0 min), 88% (3.0–4.2 min), 100% (4.21–4.5 min), 50% (4.51–5.0 min)
Flow Rate	3.0 mL/min
Column	SUPELCOSIL LC-PAH (50 mm × 4.6 mm I.D., 3 µm)
Temperature	40°C
Detection	Excitation and emission wavelength switching at each interval (Table 11)

Table 11 Wavelength Settings at Each Interval

А	Ex.: 270 nm, Em.: 330 nm
В	Ex.: 250 nm, Em.: 370 nm
С	Ex.: 330 nm, Em.: 430 nm
D	Ex.: 270 nm, Em.: 390 nm
E	Ex.: 290 nm, Em.: 430 nm
F	Ex.: 370 nm, Em.: 460 nm
G	Ex.: 270 nm, Em.: 330 nm

5-4. Amino Acid Analysis

This is an amino acid automatic analysis system with which we conducted post-column fluorescence derivatization detection using o-phthalaldehyde (OPA) / *N*-acetylcysteine as the reaction reagent. Since we used Shimadzu proprietary *N*-acetylcysteine (odorless substance) reagent, not only is handling easier than with the normally used mercaptoethanol, sensitivity for imino acids such as proline is improved. Also, because the mobile phase and reaction reagents are included in a kit, time-consuming preparation of reagents is unnecessary. Fig. 16 shows a comparison of the chromatograms and *S/N* ratios for each analyte using the former model RF-10AxL and the RF-20Axs. The results demonstrate a 4 to 5 times increase in sensitivity with the RF-20Axs.





Fig. 16 Comparison of Amino Acid Analysis Using the RF-10AxL and RF-20Axs

5-5. Repeatability in Carbamate Analysis

This is an N-methyl carbamate pesticide automatic analysis system which we used for post-column fluorescence derivatization detection with OPA as the reaction reagent. This system provides excellent accuracy in analysis of N-methyl carbamate pesticides in foods and tap water. N-methyl carbamate pesticides are often used in agricultural insecticides. Accordingly, residue standards and target values have been established in relation to various restrictions, including the Japanese Positive List System and water guality management target substances, etc. When conducting analysis by HPLC, the post-column fluorescence derivatization method is used. Shimadzu Carbamate Analysis Application System is designed specifically for analyses that comply with the individual test methods specified in the Positive List System and water quality management target substance test methods.

Fig. 17 shows the results of continuous analysis of various carbamates. From the repeatability shown for each substance in Table 16, it is clear that excellent repeatability can be obtained with this system.



Fig. 17 Results of Continuous Analysis of Carbamates

Table 14	Analytical	Conditions for	r Carbamate Analysis
----------	------------	----------------	----------------------

Mobile Phase	A: Water, B: methanol
Flow Rate	1.0 mL/min, high-pressure gradient analysis
Column	Shim-pack FC-ODS (75 mm × 4.6 mm I.D., 3 µm)
Temperature	50°C
Detection	Ex.: 339 nm, Em.: 445 nm
Injection Vol.	10 µL

Table 15 Carbamate Post-Coulmn Reaction Conditions

Reaction Reagent	A: Aqueous sodium hydroxide B: OPA, borate buffer solution containing β-mercaptopropionic acid
Flow Rate	Each 0.5 mL/min
Reaction Temp.	100°C and 50°C

Table 16 Repeatability for Each Compound (n = 6)

Compound	Area Value Repeatability (%RSD)	Compound	Area Value Repeatability (%RSD)
Aldicarb sulfoxide	2.63	Bentiocarb	2.44
Aldicarb sulfone	2.53	Carbaryl	1.62
Oxamyl	2.57	Ethiofencarb	2.69
Aldicarb	1.42	Fenobucarb	1.96

5-6. Analysis of Reducing Sugars

Various techniques are used for the separation and detection of sugars. For detection, differential refractive index detectors and evaporative light scattering detectors are commonly used. However, quantitation can be difficult with these detectors from the standpoints of separation conditions, sensitivity, and selectivity. Shimadzu Nexera Reducing Sugar Analysis System is specially designed for highly selective detection of reducing sugars. After separation of sugars by a ligand exchange, borate complex anion exchange, use of HILIC column, etc., the sugars are combined with a reaction reagent containing borate and arginine under high temperature to form fluorescent derivatives. This system allows analysis of reducing sugars in samples containing many contaminants, with excellent selectivity and sensitivity. Here we show the chromatogram of hydrolyzed cellulose nano fiber (CNF) using this system (Fig. 18).



Fig. 18 Chromatogram of Constituent Sugar Obtained from CNF

>>>> Learn More

Constituent Sugar Analysis of Cellulose Nanofibers

Nexera, Shim-pack and Prominence are trademarks of Shimadzu Corporation or its affiliated companies in Japan and/or other countries

🕀 SHIMADZU

For Research Use Only. Not for use in diagnostic procedures.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu.

Company names, products/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation, its subsidiaries or its affiliates, whether or not they are used with trademark symbol "TM" or "®". Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®".

Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.

First Edition: October, 2010 (Second Edition: August 2022)

Shimadzu Corporation www.shimadzu.com/an/