

Application News

No. L544

High Performance Liquid Chromatography

Simple and Quick Analysis of Theanine in Tea by Automatic Pre-Column Derivatization Method

Theanine, which has an *umami* taste, accounts for more than half of the free amino acids in tea and is an important component that determines the flavor of tea. It is known that there is a strong correlation between theanine content and tea quality, especially for green tea. Theanine also has various physiological effects, such as relaxing effect, stress-reducing effect, premenstrual syndrome (PMS) alleviation, and so on. With the recent health food boom, theanine is getting attention day by day.

Application News L529A introduced amino acid analysis by the automatic pre-column derivatization method using the co-injection mode of the i-Series. This article introduces an example of simple and quick analysis of theanine and other major amino acids, by the same method as L529A using the Nexera™ Series.

T. Yoshioka

Simultaneous Analysis of Amino Acids

In the Nexera Series, an automatic pretreatment function is provided in the autosampler, and its co-injection mode was used here. Primary amino acids were derivatized into fluorescent substances in the needle by using *o*-phthalaldehyde (OPA) to analyze. Fig. 1 shows the setting of the co-injection mode.

Fig. 2 shows the chromatogram of an amino acid mixed standard solution that contains theanine and 19 proteinogenic amino acids except proline. All 20 amino acids were separated within 15 min. The analytical conditions are shown in Table 1 to Table 4.




Fig. 1 Pretreatment (Co-Injection) Setting Screen of Nexera Series

Table 1 Mobile Phase Preparation Method

- 20 mmol/L (sodium) acetate buffer (pH6)
Add 2.67 g of sodium acetate trihydrate and 41 µL of acetic acid into 1 L of ultrapure water.
- 20 mmol/L (sodium) acetate buffer (pH5) containing 0.5 mmol/L EDTA-2Na
Add 0.19 g of EDTA-2Na, 2.03 g of sodium acetate trihydrate, and 308 µL of acetic acid into 1 L of ultrapure water.

Table 2 Derivatizing Reagent Preparation Method

- 0.1 mol/L borate buffer (pH9)
Add 0.62 g of boric acid and 0.20 g of sodium hydroxide into 100 mL of ultrapure water.
- MPA Reagent
Add 10 µL of 3-mercaptopropionic acid into 10 mL of 0.1 mol/L borate buffer.
- OPA Reagent
Add 0.3 mL of ethanol into 10 mg of *o*-phthalaldehyde and dissolve completely. Then add 0.7 mL of 0.1 mol/L borate buffer and 4 mL of ultrapure water.
- MPA / OPA Solution
Mix equal volume of MPA reagent and OPA reagent.

Table 3 Analytical Conditions

| | |
|------------------|--|
| Column | : Shim-pack™ XR-ODS II 100 mmL.×3 mmI.D., 2.2 µm |
| Mode | : Low pressure gradient |
| Mobile phase | : A) 20 mmol/L (sodium) acetate buffer (pH6) B) Water / Acetonitrile = 1 / 9 C) 20 mmol/L (sodium) acetate buffer (pH5) containing 0.5 mmol/L EDTA-2Na |
| Flow rate | : 1.0 mL/min |
| Column temp. | : 40 °C |
| Injection volume | : 1 µL |
| Detection | : Fluorescence detector Ex. 350 nm, Em. 450 nm |

Table 4 Time Program

| Time (min) | A.conc | B.conc | C.conc |
|------------|--------|--------|--------|
| 0 | 95 | 5 | 0 |
| 0.2 | 93 | 7 | 0 |
| 1 | 93 | 7 | 0 |
| 4 | 87 | 13 | 0 |
| 5 | 0 | 15 | 85 |
| 7.5 | 0 | 30 | 70 |
| 12 | 0 | 35 | 65 |
| 14 | 0 | 45 | 55 |
| 14.01 | 0 | 95 | 5 |
| 17 | 0 | 95 | 5 |
| 17.01 | 95 | 5 | 0 |
| 19.5 | 95 | 5 | 0 |




Fig. 2 Chromatogram of Amino Acid Mixed Standard Solution (25 µmol/L Each)

Linearity and Repeatability

Calibration curves for theanine, aspartic acid, glutamic acid, serine, glutamine, and arginine, which are the major amino acids in tea, were prepared from the chromatogram of the amino acid mixed standard solution. Calibration curves for each of these compounds were prepared in the concentration ranges of 0.25, 1.25, 2.5, 12.5 and 25 µmol/L, and their linearity (R^2) was evaluated. Area repeatability was also evaluated by a repeated analysis ($n=6$) for 12.5 µmol/L. Table 5 shows these results.

Table 5 Linearity of Calibration Curves and Repeatability of Each Amino Acid

| Amino acids | Linearity (R^2) | Repeatability (%RSD) |
|---------------|---------------------|----------------------|
| Theanine | 0.9999 | 0.37 |
| Aspartic acid | 0.9998 | 0.61 |
| Glutamic acid | 0.9997 | 0.55 |
| Serine | 0.9998 | 1.26 |
| Glutamine | 0.9999 | 0.59 |
| Arginine | 0.9999 | 0.72 |

Example of Analysis of Teas

Each of the tea-leaves was extracted by the general methods to prepare analytical samples except green tea (PET bottle). These samples were diluted 50 times with ultrapure water, filtered with a 0.22 µm membrane filter, and then used in the analysis.

Fig. 3 shows the quantitative analysis results for the amino acids in the 5 types of tea, and Fig. 4 shows the chromatograms of teas. Judging from these data, the content of amino acids varies in accordance with the type of tea.




Fig. 3 Contents of Amino Acids in Teas




Fig. 4 Chromatograms of Teas
a: Green Tea (*kabusecha*), b: Green Tea (*sencha*), c: Green Tea (PET bottle), d: Oolong Tea, e: Black Tea

Conclusion

In this analysis, derivatization of amino acids was performed in the needle by using the co-injection mode installed in the Nexera Series as a standard feature. As a result, troublesome pretreatment operation and a complex equipment configuration are not necessary. Thus, it has become possible to perform simple and quick analyses of amino acids.

Nexera and Shim-pack are trademarks of Shimadzu Corporation in Japan and/or other countries.

First Edition: Jul. 2019



For Research Use Only. Not for use in diagnostic procedure.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. Shimadzu disclaims any proprietary interest in trademarks and trade names used in this publication other than its own. See <http://www.shimadzu.com/about/trademarks/index.html> for details.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.

Shimadzu Corporation

www.shimadzu.com/an/