

Sialic Acid Stabilizing Kit for Linkage Isomer Discrimination

SialoCapper-ID Kit



SialoCapper™-ID Kit reliably stabilizes sialic acid residues, enabling linkage isomer discrimination.

Our technology prevents the decomposition of sialic acid residues and ionization suppression. Furthermore, makes it possible to discriminate the linkage isomers of sialic acids by mass spectrometry.

Features




- Simple analysis of sialic acid linkage isomer
- Improving sialoglycan sensitivity
- Versatility and scalability

Sialic Acid Linkage-Specific Alkylamidation (SALSA method)



This is our patented technique that allows stabilization by chemical modification while changing the mass of the sialic acid depending on the linkage types. As a result, the linkage isomers can be discriminated by mass spectrometry. The SialoCapper-ID Kit is a novel reagent kit for glycan pretreatment that simplifies the SALSA procedure.

Reaction scheme of SALSA method



MALDI mass spectra of labeled glycan standards derivatized with the SialoCapper-ID Kit.


Various types of mass spectrometers can be used for measuring glycans derivatized with the SialoCapper-ID Kit



The SALSA method involved in the SialoCapper-ID Kit was developed by Shimadzu corporation and improved in cooperation with Prof. Jun-Ichi Furukawa and Dr. Hisatoshi Hanamatsu at Hokkaido University.


Simple analysis of sialic acid linkage isomer

The linkage types of sialic acid can be determined from MS¹ analysis by using SialoCapper-ID Kit. This eliminates the need for separation by LC, preparation of glycan standards, and MS/MS analysis. This kit enables highly reliable measurement data acquisition and labor-saving experiments.




Improving sialoglycan sensitivity

SialoCapper-ID Kit improves sialoglycan sensitivity while stabilizing sialic acids and reducing their degradation. It also makes it easier to measure trace samples such as biological samples. The mass spectrum becomes simpler and glycan analysis becomes easier.



Versatility and scalability

SialoCapper-ID Kit can be used in combination with previous experimental procedures. Liquid-phase reaction is effective for both labeled and non-labeled glycans. In the solid-phase reaction used in combination with glycan binding beads such as BlotGlyco®, glycan purification and sialic acid derivatization can be achieved simultaneously.



Specifications : SialoCapper-ID Kit

Kit Contents	Reagent A (1 bottle), B (10 tubes), C (1 bottle)
Number of samples	Sufficient for 10 experiments/kit <ul style="list-style-type: none"> • Solid-phase reaction: Sufficient for 10 separate samples/experiment, 100 samples/kit • Liquid-phase reaction: sufficient for 50 separate samples/experiment, 500 samples/kit
Derivatization Targets	Released N-linked glycans (including labeled forms)
Storage Temperature	2-8 °C


- For solid-phase reaction, glycan purification beads are required but not included in SialoCapper-ID Kit.
Recommended: Sumitomo Bakelite BlotGlyco® BS-45414 (for 10 samples), BS-45415 (for 50 samples), BS-45407 (for 100 samples)
- For liquid-phase reaction, HILIC mode glycan purification columns or tips are required but not included in SialoCapper-ID Kit.

Glycan composition estimation software

Supporting Tool for SialoCapper-ID Kit



The Supporting Tool for SialoCapper-ID Kit can estimate glycan composition from given m/z values in the results of mass spectrometry and deal with the mass change of sialic acid residues by sialic acid derivatization such as SALSA method.



Supporting Tool for SialoCapper-ID Kit

Mass(m/z, obs.)	Intensity	Charge	Composition	Ion Species	Salt Formation	m/z(alc.)	Error(ppm)
2428.6800	0.0000		[AA-Hex3Man5GlcNAc6-1HexGlcNAc3-1]	4+		2428.6128	13.1
2456.9400 (2 matches)							
2456.9400	0.0000		[AA-Hex3Man5GlcNAc6-1GlcNAc2-2]	4+		2456.9841	1.7
2574.9700 (1 match)							
2574.9700	0.0000		[AA-Hex3Man5GlcNAc6-1HexGlcNAc3-1]	4+		2574.9707	0.3
2603.0000 (1 match)							
2603.0000	0.0000		[AA-Hex3Man5GlcNAc6-1GlcNAc2-2]	4+		2603.0009	0.8
2748.9900 (1 match)							
2748.9900	0.0000		[AA-Hex3Man5GlcNAc6-1HexGlcNAc3-1]	4+		2748.9347	16.2
2777.0200 (1 match)							
2777.0200	0.0000		[AA-Hex3Man5GlcNAc6-1HexGlcNAc3-1]	4+		2777.0660	13.0
2822.0400 (2 matches)							
2822.0400	0.0000		[AA-Hex3Man5GlcNAc6-1HexGlcNAc3-1]	4+		2822.0510	3.6
2822.0400	0.0000		[AA-Hex3Man5GlcNAc6-1GlcNAc2-2]	4+		2822.0763	12.9
3114.0700 (1 match)							
3114.0700	0.0000		[AA-Hex3Man5GlcNAc6-1HexGlcNAc3-1]	4+		3114.0609	33.3
3435.1400 (1 match)							
3435.1400	0.0000		[AA-Hex3Man5GlcNAc6-1HexGlcNAc3-1]	4+		3435.2729	38.7

Glycan composition candidates

Recommended specifications for computers

- Operating system: Windows® 10 Pro 64 bit (Japanese/English)
- Memory: More than 8 GB
- Storage: More than 256 GB
- Optical drive: DVD drive

*A computer is not included for installing the software.

SialoCapper is a trademark of Shimadzu Corporation.
BlotGlyco is a registered trademark of Sumitomo Bakelite Co., Ltd.
Windows is either a registered trademark or a trademark of Microsoft Corporation in the United States and/or other countries.



Shimadzu Corporation

www.shimadzu.com/an/

For Research Use Only. Not for use in diagnostic procedures.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.
Company names, products/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation, its subsidiaries or its affiliates, whether or not they are used with trademark symbol "TM" or "®".
Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®".
Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The contents of this publication are provided to you "as is" without warranty of any kind, and are subject to change without notice. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication.