# Isolate Cannabinoids using





• THC and CBD interact to different degrees with the C18 Media

• As a result CBD passes through the media faster and exits the column first

### THC is collected later





### **Gallons of solvent to**

# **Biotage systems**

	Capacity/24hr	Gallons of solvent used/24hr
Isolera <sup>TM</sup> LS	1.2 kg	50 - 70
Biotage® Flash 150	3.6 kg	150 - 200
Biotage® Flash 400	28.8 kg	1200 - 1500



# produce purified distillate \*\*





Your results will vary



**F** 

\$ 1.04

			\$ 0.90	
Cost per Gram				
	50 grams per run Isolera™LS	250 grams per run Biotage® Flash 150	2000 grams per run Biotage® Flash 400	
	* Cost per rui	n assumes that	t columns are	

re-used 50 times before replacement

## Want to go more in-depth?



QR Code to THC Removal Blog



QR Code to White Paper on Flash

